

DRIPASSIST™

Infusion Rate Monitor



Simple saves lives.



DRIPASSIST™

Infusion Rate Monitor

Redefining the standard of care for gravity infusion therapy.

SAFE, ACCURATE IV
MEDICATION MANAGEMENT

The award-winning DripAssist Infusion Rate Monitor provides easy, reliable IV rate monitoring at an affordable price.

DripAssist quickly attaches to any gravity drip set to continuously watch flow rate and volume for fluids and medication.



💧 Easy

DripAssist does the drop counting and calculating so you can focus on your patient. Takes less than five minutes to learn.

💧 Reliable

Get fast, precise measurement for your IV infusions. Measure drip rate with 99% accuracy. Runs continuously for two weeks on one AA battery.

💧 Affordable

Monitor more patients with less expense. Costs a fraction of the price of an infusion pump and works with any drip set, with no proprietary consumables to buy.

- US FDA cleared
- Assembled in USA
- No calibration or maintenance
- Works with most IV drip sets
- Uses one AA battery
- Learn in less than 5 minutes





Safer Infusions in Any Care Environment

Clinicians across the globe use DripAssist for safe, accurate IV medication dosing in a wide range of settings.

Outpatient and Ambulatory Infusion - Saving Nursing Time and Reducing Costs

CHEMOTHERAPY, ANTIBIOTICS, IVIG, CHRONIC
INFLAMMATORY DISEASE

Maternal Care - Safer Labor and Delivery

MAGNESIUM SULFATE, OXYTOCIN

Pediatric and Elderly Care

FLUID MANAGEMENT

Pre-Hospital Medication Administration

EPINEPHRINE, VASOPRESSORS, NITROGLYCERIN

Infectious Disease Management

QUININE, FLUID MANAGEMENT FOR CHOLERA

CLINICAL TESTIMONIALS

"I used DripAssist with a 16 month old oncology patient receiving platelets and I just loved it. Easy to setup and use. It was very accurate and easy to control the drip rate to what I wanted. This will be put to great use."

PEDATRIC INFUSION COORDINATOR,
NEW HANOVER MEDICAL CENTER, USA

"As a center where we do deliveries and don't have a large staff, DripAssist ensures that the patients are getting the correct dose of solutions. The device acts as another staff that will never sleep and always be present."

DR. EDOUARD BORDES, MEDICAL DIRECTOR
COMMUNITY CLINIC IN PORT-AU-PRICE, HAITI

"DripAssist is a novel product for agencies seeking a solution that decreases setup time, improves accuracy and is less expensive than a pump."

SAMUEL KORDICK, CLINICAL SPECIALIST,
CYPRESS CREEK EMS, USA

AWARD WINNING DESIGN SUPPORTED BY USAID AND USED IN OVER 15 COUNTRIES



**2016 Named One of 12 Most Important Health
Innovations by Popular Science Magazine**



**2015 International Design Excellence Awards
from Industrial Design Society of America**



Learn how to get started with DripAssist. Visit www.shiftlabs.com

DRIPASSIST™

Infusion Rate Monitor



Specifications

Battery	One AA battery, lasts 360 hours
Weight	108 g (3.8 oz)
Size	127 x 62 x 28mm (5 x 2.4 x 1.1in)
Measurement Units	Flow rate (mL/h), drops per minute (dp/m) and total volume (mL)
Alarm	80 dB at 10cm approx. Sounds at $\pm 13\%$ rate change or when flow stops
Minimum Chamber Diameter	13 mm (0.51in)
Maximum Chamber Diameter	21.1 mm (0.83in)
Compatible Drip Sets	10, 15, 20, 60 gtt/mL
Accuracy	$\pm 1\%$ drip rate
Warranty	1 year limited warranty

Standards Compliance

Radiated Emissions	CISPR 11:2010
Ingress Protection	IP22
Electromagnetic compatibility	IEC 61000-4
Regulatory	US FDA 510(k) K150687, CE Mark, ISO 13485:2003 and ISO 13485:2003 CMDCAS

Additional Warnings for Customers in the USA

The following warnings apply only for customers intending to use the DripAssist device within the United States of America:

- Only trained, qualified personnel should operate the DripAssist.
- The DripAssist is not intended for home use.
- Caution: Federal law restricts this device to sale by or on the order of a licensed practitioner.



www.shiftlabs.com
info@shiftlabs.com
+1 206.369.1216

Corporate Headquarters
1600 Dexter Ave N • Suite A
Seattle, WA 98109
United States

Mailing Address
1752 NW Market St. #211
Seattle, WA 98107
United States

DRIPASSIST™

Infusion Rate Monitor

WELCOME

Getting Started with DripAssist

DripAssist is a simple, accurate Infusion Rate Monitor that easily clips onto any gravity drip set to **continuously monitor** flow rate and volume when administering IV medications.

How is DripAssist different from a pump?

- DripAssist does not control an IV infusion. It is a companion technology for manual, gravity infusion.
- To use DripAssist, the clinician controls the flow with the roller clamp attached to a tubing set.
- DripAssist monitors that flow and, when the alarm is activated, can alert the user to changes in the set flow rate.

For what types of care should DripAssist be used?

DripAssist can be used with any gravity infusion and using any standard tubing set with a gtt/mL factor of 10, 15, 20, or 60 gtt/mL. DripAssist is not a replacement for a pump.

Where can I find best practice instructions for DripAssist?

We offer a series of videos to assist you in operating your DripAssist. You can also download an app for Android called **Infusion Coach** which includes information on how to use DripAssist, best practices for gravity infusion, and detailed information on cleaning, positioning, and using the alarm.

How long does the battery last with continual use?

If you use DripAssist for 24 hours a day, the AA battery will last approximately two weeks, or 290 hours. You may experience some variation based on the quality of the batteries used.

How do I clean my DripAssist?

Clean your DripAssist using a standard cleaner or disinfectant. Follow the instructions on the container. For more information, see our video “How To Clean DripAssist”.